



Project Management and Software Development
for Medical Applications

Image Processing for Digital Breast Tomosynthesis

Intermediate Presentation, 20.06.2023

Project Supervisor: Mickael Tardy

Presented by: Berfin Kavsut



Technische Universität München



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING



Section 01

Background

Medical Imaging Modality

Digital Breast Tomosynthesis [1]

- 3D images
- Reconstructed from 2D X-ray projections
- Cross-sectional visualization via slices




Figure 1. Acquisition Geometry of Digital Breast Tomosynthesis [2]

[1] "Breast Imaging Diagnosis and Screening in Korea," *Breast Imaging, 12th International Workshop IWDM*, Gifu City, Japan, 2014, pp. xxvii-xxviii, doi: 10.1007/978-3-319-07887-8.

[2] D. Kontos, P. R. Bakic, and A. D. A. Maidment, "Texture in digital breast tomosynthesis: a comparison between mammographic and n tomographic characterization of parenchymal properties," in *Proc. SPIE*, vol. 6915, Mar. 17, 2008, p. 69150A. doi: 10.1117/12.773144.



Slabbing


Slices:

- 1 mm thickness for standard protocol
- very high image noise level [3]

Slabs:

- 10 mm thickness for standard protocol
- better image quality [3]

"thin" slices
separated by 1 mm



fewer "thick" over-
lapping slices



Figure 2. Slabbing slices [4]

[3] F. Diekmann, H. Meyer, S. Diekmann et al., "Thick Slices from Tomosynthesis Data Sets: Phantom Study for the Evaluation of Different Algorithms," *Journal of Digital Imaging*, vol. 22, no. 5, pp. 519-526, Oct. 2009, doi: 10.1007/s10278-007-9075-y.

[4] HealthManagement.org, "Digital Breast Tomosynthesis in screening – approaches to reduce reading time," HealthManagement, <https://healthmanagement.org/c/decision-support/whitepaper/digital-breast-tomosynthesis-in-screening-approaches-to-reduce-reading-time> (accessed May 13, 2023).





Section 02

Methods



Technische Universität München



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING

Pipeline




Figure 3. Pipeline schematic of the project

Pre-Processing

- Background noise removal
- Clearing the bright breast skin




Figure 4. Breast skin segmentation

Pre-Processing

Breast segmentation

Steps:

1. Breast segmentation with **triangle thresholding**
2. **Lower resolution** of the segmented breast mask
3. Segment breast skin by using the **morphological operation** called erosion
4. Create breast, skin, and background masks




Figure 5. Triangle Thresholding [5]



[5] "Image thresholding," OpenCV, https://docs.opencv.org/3.4/d7/d4d/tutorial_py_thresholding.html (accessed Jun. 19, 2023).

Projection Methods

- Maximum intensity projection: $f_w^{MIP}(x) = 1 \text{ if } x = 1, \text{ else } 0$
- Average intensity projection: $f_w^{AIP}(x) = 1$
- SoftMIP: $f_w^{softMIP}(x) = x^4$

create projection line: $P_s(x), 0 \leq x \leq l$
 $P_s(0) = \min(P), P_s(l) = \max(P).$

compute projection: $p = \frac{1}{\int_0^1 f_w(x) dx} \int_0^1 f_w(x) P_s(x) dx$




Figure 6. Maximum intensity projection and average intensity projection [6]



Section 03

Software Development



Technische Universität München



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING

Software Tools

Version control: GitHub

IDE: PyCharm

Issue Tracking System: GitHub issues



Python Libraries:

- scikit-image
- OpenCV
- pandas
- pydicom
- numpy



Issue Tracking

To do | 7

Keep track of DICOM's
SeriesInstanceUID when saving slabs

Save config as json alongside with
generated slabs

Provide tools allowing for processing
a DICOM DBT volume as input and
saving slabs as output

Choose the configuration method to
call the slab generator

Displaying pathology better (related
to normalization methods)

Normalization when creating the
slabs

...

In progress | 4

Do not apply averaging on the skin
border pixel when skin-border is
activated

Add configuration parameter to slab
generator to skip slices near to ends

Provide default config values to
avoid failures

Challenge breast skin effect by
choosing dicom files with thick and
bright breast skins

+

Done | 6

Read box coordinates of pathologies
(20 samples)

Use lower resolution images for
breast segmentation

Add SoftMIP projection

Saved 10-bit images with 16-bits
range

Allow DICOMFileProcessor to use
images from local storage

Choose the morphological structure

+





Section 04

Challenges



Technische Universität München



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING

Challenges

Design the end-to-end architecture

Run code only once

- Prepare dataset by downloading DICOM files
- Generate multiple configurations
- Save the results with configuration

Projection with and w/o breast skin

- Keep the breast skin from center image
- Modify projection accordingly





Section 05

Gantt Chart




Technische Universität München




JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING

Gantt Chart



Updated Gantt Chart





Section 06

Appendix



Technische Universität München



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING



UML Class Diagram

- Current design



Technische Universität München



JOHNS HOPKINS

WHITING SCHOOL
of ENGINEERING

